## Mouse Anti-HCV Core Protein Monoclonal Antibody Datasheet

Product Name: mAb anti-HCV Core Protein Clone No.: A1/3D1

Catalogue No.: MO-I40015A Quantity: 0.5 mg/vial

**Description:** Mouse monoclonal antibody to human

hepatitis C virus (HCV) core protein

Purification: Protein G affinity purified

**Product Type:** Primary antibody

Target Protein: Human hepatitis C virus (HCV) core protein

Immunogen: Synthetic peptides derived from HCV core

protein

Fusion Sp2/0-Ag14

Myeloma:

**Specificity:** mAb A1/3D1 is reactive to recombinant core

protein C + envelope protein M (residues 1-142 on HCV polyprotein) and synthetic core

protein C (residues 1-61 on HCV

polyprotein).

**Species** Human hepatitis C virus, others not tested

Reactivity:

**Cross -** It showed no cross reaction with

**Reactivity:** recombinant or synthetic HCV non-structural

proteins (NS-3 and NS-4).

Host / Isotype: Mouse, IgG1 Kappa

Formulation: Lyophilized from a solution in 0.01M PBS, pH

7.2

**Reconstitution:** Double distillated water is recommended to

adjust the final concentration to 1.00mg/mL.

Storage: Store at -20°C

Research Virology

Area:

**Background:** Hepatitis C virus (HCV) causes chronic

hepatitis and liver cirrhosis in human through blood and body fluid transmission. HCV has a positive sense single RNA genome enclosed in the nucleocapsid made of Core Protein (Capsid Protein). The nucleocapsid

is covered by an envelope made of lipoproteins (E1 and E2). The 9.6 kb HCV genome has a single open-reading frame, which is to be translated into a single polyprotein. HCV viral proteins are

produced after processing the polyprotein.

Genes for core protein and envelop proteins are located adjacently at the 5'-end of HCV genome, followed by genes for non-structural proteins including NS2, NS3,

NS4A, NS4B, NS5, NS5A and NS5B.

**Applications:** ELISA: mAb A1/3D1 was tested using

indirect ELISA method. The testing plates were coated with HCV non-structural protein (NS), core protein (C) and mixture of NS and C respectively. The neat culture supernatant of hybridoma showed a strong reactivity with HCV capsid protein (OD492 value > 2.0) and with mixture of NS and C (OD492 value >1.5), and showed no reactivity with HCV non-structural protein (OD492 value < 0.03). mAb A1/3D1 recognizes in-vitro translated HCV core

protein

Not suitable for use in Western Blot

**References:** If research is published using this product,

please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same

category as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture ONLY, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.