ANOGEN - A Division of YES Biotech Laboratories Ltd. 2355 Derry Road East, Unit 23, Mississauga, ON, Canada L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

## Mouse Anti-Human MIP-1 $\alpha$ Monoclonal Antibody Datasheet

Product Name: mAb anti-human MIP-1α Catalogue No.: MO-C40051N		Clone No.: 376 Quantity: 0.5 mg/vial	
Purification:	Protein G affinity purified	Background:	Macrophage inflammatory protein 1α (MIP-1α), also known as CCR3, belongs
Product Type:	Primary antibody		to the MIP-1 CC chemokine family.
Target Protein:	Human MIP-1α		The cytokine is mainly produced by macrophages after stimulation by
Immunogen:	Purified recombinant human MIP-1 $lpha$		bacterial endotoxins. In acute inflammatory stage, MIP-1α plays an
Fusion Myeloma:	Sp2/0-Ag14		important role in the recruitment and activation of granulocytes. MIP-1 $\alpha$
Specificity:	Reactive with recombinant human MIP-1 $\alpha$		also stimulates the production of pro- inflammatory cytokines, including IL-1, IL-6 and TNF-α. In normal status, MIP- 1α plays a role in promoting
Species Reactivity:	Human, others not tested		homoeostasis.
Cross-	Does not show any cross reaction with	Applications:	ELISA, Western blot
Reactivity:	other human cytokines or growth factors tested such as M-CSF, GM-CSF, IL-7, IL-8, IL-16, IL-1 $\beta$ , bFGF, TGF- $\beta$ , EGF, TNF- $\alpha$ and EPO.	References:	If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one
Host / Isotype:	Mouse, IgG1 Kappa		unit of product in the same category as gratitude.
Formulation:	Lyophilized from a solution in 0.01M PBS, pH 7.2		
Reconstitution:	Double distilled water is recommended to adjust the final concentration to 1.00mg/mL.		

This product is for LAB	ORATORY RESEARCH USE and further manufacture ONLY, and cannot be administrate	d to human and
	animals for use in diagnostic and therapeutic procedures.	
	Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.	
Page 1 of 1 page(s)		S7.5 (02)
	To place order, please contact us by phone, fax or by email: <u>info@anogen.ca</u> ,	
	or with our secure online store: <u>www.anogen.net</u>	