



Mouse Anti-Human IL-17A Monoclonal Antibody Datasheet

Product Name: mAb anti-Human IL-17A

Clone No.: 2F3

Lot No.:

Catalogue No.: MO-C40093B

Quantity: 0.5 mg

Expiration:

Description: Mouse monoclonal antibody to human interleukin-17A (IL-17A)

Product Type: Primary antibody.

Purification: Protein G affinity purified

Target Protein: Human IL-17A

Immunogen: Purified recombinant human IL-17A

Fusion Myeloma: SP2/0

Specificity: Reactive with human IL-17A

Species Reactivity: Human, others not tested

Host / Isotype: Mouse, IgG2b, Kappa

Formulation: Lyophilized from a solution in 0.01M PBS, pH 7.0

Reconstitution: Double distilled water is recommended to reconstitute the antibody to 1mg/ml

Storage: Store at -20°C

Research Area: Cytokine. Inflammatory and auto-immune diseases


Background: IL-17A is an inflammatory cytokine produced by Th17 cells, a newly characterised group of CD4 T help cells that plays a role in autoimmune diseases and chronic inflammatory diseases. IL-23 is the main inducer of Th17 cell development. Other cytokines such as TGF-beta, IL-6 and IL-21 were also found to contribute to the procedure. IL-17A belongs to the IL-17 family which includes IL-17A, IL-17B, IL-17C, IL-17D and IL-17E. IL-17, like interferon-gamma, is a

potent mediator of delayed-type reaction. IL-17 induces the secretion of many other pro-inflammatory cytokines and chemokines such as IL-6, IL-8, GM-CSF, G-CSF, TGF-beta, TNF-alpha, GRO-alpha, and MCP-1. IL-17 stimulates the expression of NF-kappaB, mitogen activated protein kinase, and prostaglandins by many cell types, and enhances the production of nitric oxide.

Applications:

ELISA: The anti-human IL-17A antibody clone 2F3 reacts with human IL-17A in ELISA.

Western blot: The image below shows that anti-human IL-17A antibody 2F3 reacts with human 293 cell expressed human IL-17A antigen (10ng, 8ng) in Western Blot. The blot was stained with HRP substrate 4-chlorol-1-Naphthol.



References:

If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same category as gratitude.

This product is for LABORATORY RESEARCH USE and further manufacture ONLY, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.