ANOGEN - A Division of YES Biotech Laboratories Ltd. 2355 Derry Road East, Unit 23, Mississauga, ON, Canada L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

## Mouse Anti-Human IL-17A Monoclonal Antibody Datasheet

		Clone No.: 8H7	Lot No.: Expiration:
		uantity: 0.5 mg	
Description:	Mouse monoclonal antibody to human interleukin-17A (IL-17A)		potent mediator of delayed-type reaction. IL-17 induces the secretion of many other
Product Type:	Primary antibody.		pro-inflammatory cytokines and chemokines such as IL-6, IL-8, GM-CSF, G-CSF, TGF-beta,
Purification:	Protein G affinity purified		TNF-alpha, GRO-alpha, MCP-1. IL-17 stimulates the expression of NF-kappaB,
Target Protein:	Human IL-17A		mitogen activated protein kinase, and prostaglandings by many cell types, and
Immunogen:	Purified recombinant human IL-17A		enhances the production of nitric oxide.
Fusion Myeloma:	SP2/0	Applications:	<b>ELISA:</b> The anti-human IL-17A antibody clone 8H7 reacts with human IL-17A in ELISA.
Specificity:	Reactive with human IL-17A		Western blot: The image below shows that
Species Reactivity:	Human, others not tested		this anti-human IL-17A antibody 8H7 reacts with human 293 cell expressed human IL-17A antigen (250ng) in Western Blot. The blot
Host / Isotype:	Mouse, IgG1, Kappa		was stained with HRP substrate 4-Chloro-1- naphthol.
Formulation:	Lyophilized from a solution in 0.01M PBS, p 7.0	н	MWM rlL-17A 245 KD - 180 KD -
Reconstitution:	Double distilled water is recommended to reconstitute the antibody to 1mg/ml		135 KD - 100 KD - 75 KD -
Storage:	Store at -20°C		63 KD -
Research Area:	Cytokine. Inflammatory and auto-immune diseases		48 KD - 35 KD - 25 KD - 20 KD -
Background:	IL-17A is an inflammatory cytokine produce by Th17 cells, a newly characterised group o		17 KD - 11 KD -
	CD4 T help cells that plays a role in	References:	If research is published using this product,
	autoimmune diseases and chronic		please inform Anogen in order to cite the
	inflammatory diseases. IL-23 is the main		reference on this datasheet. Anogen will
	inducer of Th17 cell development. Other		provide one unit of product in the same
	cytokines such as TGF-beta, IL-6 and IL-21		category as gratitude.
	were also found to contribute to the	ihz	
	procedure. IL-17A belongs to the IL-17 fam		
	which includes IL-17A, IL-17B, IL-17C, IL-17I and IL-17E. IL-17, like interferon-gamma, is		
This product i	s for LABORATORY RESEARCH USE and furth animals for use in diagno		
Page 1 of 2 p	Manufactured by ANOGEN - A I	Division of YES Biotech	Laboratories Ltd. S7.5 (01)

To place order, please contact us by phone, fax or by email: info@anogen.ca, or with our secure online store: www.anogen.net