HRP Conjugated Mouse Anti-Human MCP-1 Monoclonal Antibody Datasheet

Product Name: HRP conjugated mAb anti-Human MCP-1 Clone No.: \$101

Catalogue No.: MO-C40021T Quantity: 0.2 mL/vial

Description: Horseradish peroxidase (HRP) conjugated

mouse monoclonal antibody to human Monocyte Chemotactic Protein-1 (MCP-1)/ Monocyte Chemotactic and Activating

Factor (MCAF)

Purification: Protein G affinity purified

Product Type: Tracer antibody in matched antibody pair.

HRP conjugated

Target Human MCP-1

Protein:

Immunogen: Purified recombinant human MCP-1

Fusion Sp2/0-Ag14

Myeloma:

Specificity: This antibody reacts with natural and

recombinant human MCP-1.

Species Human, others not tested

Reactivity:

Cross- This antibody does not react with human reactivity: interleukin-8 (IL-8) and other human

cytokines tested such as interleukin-1 β (IL-1 β), serum amyloid A (SAA) and epidermal

growth factor (EGF).

Host / Mouse, IgG1 Kappa

Isotype:

Storage Buffer Mixture of 50% glycerol and 50% PH7.2

Formulation: 0.01M PBS. Final PH = 7.0 ± 0.1 .

Storage: Store at -20°C

Research Cytokine, chemotaxis and inflammation

Area:

Background: Monocyte chemotactic and activating

factor (MCAF) is also called monocyte chemotactic protein-1 (MCP-1) and chemokine (C-C motif) ligand 2 (CCL2). It is primarily secreted by monocytes, macrophages and dendritic cells. This cytokine displays chemotactic activity for monocytes, T-cells, and basophils, but not for neutrophils or eosinophils. MCAF causes the degranulation of basophils and mast cells, and augments the activity of monocyte and macrophage. MCAF plays an important role in inflammation, angiogenesis, auto-immune diseases,

granuloma formation.

Applications: ELISA: The HRP conjugated antibody can

be used for quantitative detection of human MCP-1 in sandwich ELISA in combination with capture antibody S14

renal diseases, chronic infection and

(Cat. No.: MO-C40021B).

References: 1. Randolph, G. J. and Furie, M. B. A

soluble gradient of endogenous monocyte chemoattractant protein-1 promotes the transendothelial migration of monocytes in vitro. J. Immunol. 1995 **155**:3610-3618.

2. Gwendalyn J. Randolph and Martha B.

This product is for LABORATORY RESEARCH USE and further manufacture ONLY, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.

Page 1 of 2 page(s)

S7.5 (03)

To place order, please contact us by phone, fax or by email: info@anogen.ca, or with our online store: www.anogen.ca



Furie. Mononuclear phagocytes egress from an in vitro model of the vascular wall by migrating across endothelium in the basal to apical direction: role of intercellular adhesion molecule 1 and the CD11/CD18 integrins. J Exp Med. Feb 1, 1996; 183(2): 451–462.

3. Mehrdad Baghestanian et al. The c-kit Ligand Stem Cell Factor and Anti-IgE Promote Expression of Monocyte Chemoattractant Protein-1 in Human Lung Mast Cells. Blood. December 1, 1997 vol. 90 no. 11 4438-4449.

If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same category as gratitude.