



# ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON, Canada L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

## Human Chorionic Gonadotrophin (hCG) Datasheet

---

**Product Name:** Human Chorionic Gonadotrophin (hCG)    **Catalogue No.:** AP30002    **Quantity:** 10 mg/vial

---

**Description:** HCG from human urine

**Source:** Human Urine

**Formulation:** Lyophilized powder

**Product Form:** Lyophilized powder. Use distilled water to dilute.

**Reconstitution:** Double distilled water is recommended to adjust the final concentration to 1.00 mg/mL.

**Stability:** Should be stored below 0 °C. Reconstituted form is best stored refrigerated at 4 °C.

**Activity:** Biopotency: 3200IU/mg.  
Immunopotency: 6000IU/mg.

**References:** If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same category as gratitude.

This product is for **LABORATORY RESEARCH USE** and further manufacture **ONLY**, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.