



ANOGEN - A Division of YES Biotech Laboratories Ltd.

2355 Derry Road East, Unit 23, Mississauga, ON, Canada L5S 1V6 • Tel: (905) 677-9221 • Fax: (905) 677-0023

Human Umbilical Cord Serum Datasheet

Product Name: Human Umbilical Cord Serum	Catalogue No.: S30010	Quantity: 10 mL/vial
---	------------------------------	-----------------------------

Description: Human umbilical cord serum is prepared from new born placenta umbilical cord which is non-reactive to hepatitis B surface antigen, HIV-1 antibodies. The AFP content was > 0.05 mg/mL. Please note: For cell culture and other applications that require sterile condition, sterile filtration is required.

Source: Human blood

Formulation: Liquid

Storage: Short Term: 4 °C for a few days.
Long Term: -20 °C.

References: If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same category as gratitude.

This product is for **LABORATORY RESEARCH USE** and further manufacture **ONLY**, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.