



Mouse Anti-CA 19-9 Monoclonal Antibody Datasheet

Product Name: mAb anti- human carbohydrate antigen 19-9 (CA 19-9).

Clone No.: 1A5

Catalogue No.: MO-T40015C

Quantity: 0.5 mg/vial

Description: Mouse monoclonal antibody to human CA 19-9 antigen.

repeated freeze and thaw cycles.

Purification: Protein G affinity purified

Storage: Store at -20°C

Product Type: Primary antibody

Research Area: Oncology

Target Protein: Human CA 19-9

Background: CA19-9 has been the most valuable serum biomarker in monitoring progress and regression of pancreatic cancer and colorectal cancer. CA19-9, also known as sialyl-Lewis^A, is a tetrasaccharide which is expressed in abundance on the surface of cancer cells. CA 19-9 usually ranges from 0 to 37 units per milliliter in a healthy person. CA 19-9 levels can be higher in patients with pancreatic cancer and other cancers.

Immunogen: Human colon adenocarcinoma cell line SW1116

Fusion Myeloma: Sp2/0-Ag14

Specificity: Reactive with SW1116 and also reactive with highly purified human CA19-9 antigen from the cell culture supernatant of a human Colon adenocarcinoma.

Applications: **Indirect ELISA:** This antibody at a wide concentration range can react with ELISA plate coated with SW1116 and also with ELISA plate coated with purified CA19-9 at 0.125 KU per well.

Reactivity: Human, others not tested

Species Reactivity: Not available

Host / Isotype: Mouse, IgG2a Kappa

References: If research is published using this product, please inform Anogen in order to cite the reference on this datasheet. Anogen will provide one unit of product in the same category as gratitude.

Formulation: Lyophilized from a solution in 0.01M PBS, pH 7.2

Reconstitution: Double distilled water is recommended to reconstitute the antibody. Avoid

This product is for **LABORATORY RESEARCH USE** and further manufacture **ONLY**, and cannot be administrated to human and animals for use in diagnostic and therapeutic procedures.

Manufactured by ANOGEN - A Division of YES Biotech Laboratories Ltd.